

# Iodixanol Gradient Scale Up Using Large Scale Ultracentrifugation

Sandra Meriño

Director of Global Strategic Solutions

**KII & PKII**  
Continuous Flow Ultracentrifuges

Alfa Wassermann Promatix 1000™ Bio-purification system




Alfa Wassermann AFH-KII automated ultracentrifuge is a scalable centrifuge that is used in cGMP manufacture of gene therapy products



# Laboratory scale ultracentrifuges allow a vast range of purifications to be made at a small scale – but these are not scalable for manufacturing

These rotors are limited by their capacity; too many runs would be needed to make a manufacturing batch of product at sufficient volume

**Principle:** Fill tube – spin tube – collect from tube



Type of Rotor	Pelleting	Type of separation		
		Rate-zonal sedimentation	Rate-zonal flotation	Isopycnic
Fixed-angle	Excellent	Limited	Good	Variable*
Near vertical	NS	Poor	Good	Variable*
Vertical	NS	Good	Good	Excellent
Swinging bucket	Inefficient	Good	Excellent	Good**
Zonal	NS	Excellent	Excellent	Good

\* good for macromolecules, poor for cells and organelles

\*\*good for cells and organelles, caution needed if CsCl used

NS – not suitable


D. Rickwood, T.C. Ford and J. Steensgaard, 1994 Centrifugation Essential Data.




# Scaling up ultracentrifugation usually operates in continuous flow mode

These rotors (both disc and tubular types) are not limited by their size as fluid continuously enters and leaves the rotor during high speed operation.

Principle: Load Gradient – Spin Rotor – Flow Product – Brake Rotor – Collect Gradient





Parameter	CF32	AFH-KII
Shape	Disc type	Tubular
Speed max.	32 000 rpm	40 500rpm
Centrifugal force	102 000 xg	121 200 xg
Capacity	430 mL	3200 mL
K factor	42	29.7
Flow Path	Loop	Dual inlet
Scalable	No	Yes
Automated	Manual process	Automated




Transfer from small to large scale gradients requires that the separate tubes of a rotor are replaced by one tube of liquid which fills the rotor

Tube Rotor




228 ml




12 ml (x20)


K3 Rotor



3200 ml



DECCELERATING




AT REST

160 ml (x20)

The gradient itself reorients, within the rotor during the acceleration and deceleration phases



# Scale up is made using the rotor technology of creating different volume separations within the same rotor



Maximum radius – must not change  
 Minimum radius – must not change  
 Gradient load density – must not change

Gradient load volumes – increases volumetrically on scale up  
 Product flow rate – increases volumetrically on scale up



Particle Range	700S
Flow / Batch	Continuous Flow
Rotor Type	K3 3200 PK3-1600 PK3-800 PK3-400 PX3-230 PX3-120
Flow Rate	28 L/h 14 L/h 7 L/h 3.6 L/h 2.0 L/h 1.0 L/h



# Alfa Wassermann ultracentrifuge system capacities



Promatix

	<u>AW Promatix 1000®</u>	<u>PKII</u>	<u>KII</u>
<b>Typical Feed Flow</b>	0.25 – 2 L/h	up to 15 L/h	up to 30 L/h
<b>Rotor Size</b>	PX3 – 120 mL PX3 – 230 mL	PK3 – 400 mL PK3 – 800 mL PK3/PK6 – 1600 mL	K3/K6 – 3200 mL <i>K10 – up to 8 Liters</i> <i>K5 – 8.4 Liters</i>
<b>Batch Volume (5h run)</b>	Up to 5 L	Up to 75 L	Up to 150 L
<b>Max. Rotation Speed</b>	35,000 rpm	40,500 rpm	40,500 rpm
<b>Gravitational Forces</b>	Up to 90,500xg	Up to 121,200xg	Up to 121,200xg
<b>Scale Factor</b>	27x scale down 14x scale down	8x scale down 4x scale down 2x scale down	1x scale



PKII



KII



Scale up efficiency means that large volumes of material can be processed in a short time

Scale Volume	Process Volume	Rotor Volume	Run operations / days	Process 1 liter
Laboratory	23ml (10%)	38ml x 6 = 228ml	16 h (overnight)	44 runs
Laboratory Batch	160ml (10%)	1600ml	16 h (overnight)	6 runs
Laboratory Continuous flow	Up to 20 liters (5 L/h)*	430ml	8 h (day run)	12 min
AW PKII	Up to 25 liters (3.7 L/h)*	400ml	8 h (day run)	2 min
AW PKII	Up to 50 liters (7.5 L/h)*	800ml	8 h (day run)	2 min
AW PKII	Up to 100 liters (15 L/h)*	1600ml	8 h (day run)	2 min
AW KII	Up to 200 liters (30 L/h)*	3200ml	8 h (day run)	2 min

\*dependent on particle size





# AAV purification using continuous flow Iodixanol density gradient ultracentrifugation in the Promatix 1000

**Bac-inCapX-inRep2**

MOI = 10  
X: any serotype

**Bac-Gene of Interest**

MOI = 5



ESF921 (serum free)

Sf9


↓ 3 days

Harvest cell pellet and prepare lysate

↓  
Continuous flow ultracentrifugation

↓  
Buffer exchange and qPCR

AAV and ultracentrifugation work by  
Dr Haifeng Chen, Virovek Inc, Hayward, CA,  
USA




SEPARATION TECHNOLOGIES




# Reorienting density gradient centrifugation using continuous flow

1. With a stationary rotor the AFH loads the gradient.
2. Slow acceleration reorients the gradient.
3. At 35 000 rpm flow sample into the rotor; vector is captured in the density gradient, waste flows out.
4. Isopycnic sedimentation occurs in zones where the gradient density equals a particle's buoyant density.
5. Slow deceleration reorients the gradient and retains the separation zones.
6. AFH collection and analysis.



# AAV Purification Results



# Scale up of Iodixanol density gradients 27 fold using the AW scalable ultracentrifuges

**Promatix 1000  
R&D**



**PKII  
Pilot Scale**



**KII  
Manufacturing Scale**



Fixed Parameters:  
Maximum radius  
Minimum radius  
Gradient load density

Variable Parameters:  
Gradient load volumes –  
increases volumetrically

Product flow rate – increases  
volumetrically




# Scale up of Iodixanol density gradients 27 fold using the AW scalable ultracentrifuges

All systems: run at 35000rpm, 90 500xg, AAV residence time is 11.6 minutes, banding time is 2 hours.


Core/Vol.	25% (w.v) Iodixanol	50% (w.v) Iodixanol	Sample Volume	Flow Rate
PX-120mL	20 mL	45 mL	80 mL	5 mL/min
PX-230mL	40 mL	90 mL	160 mL	10 mL/min
PK-400mL	70 mL	165 mL	280 mL	17 mL/min
PK-800mL	140 mL	330 mL	560 mL	34 mL/min
PK-1600mL	280 mL	660 mL	1120 mL	68 mL/min
KII-3200mL	560 mL	1320 mL	2240 mL	136 mL/min



# Iodixanol density gradients using continuous flow ultracentrifugation scale up in a predictable pattern



It was established from the work with the AAV vector that the separation zone of the AAV in the density gradient was between 32% and 48% Iodixanol. Using this information analysis of the separation zone within this zone was made.



# Conclusion

---

- Density Gradient ultracentrifugation is a scalable process for use in gene therapy production
- Linear scale up of Iodixanol density gradients is possible

